



Chemical	Used On	How
<p>Oxivir® 1 Wipes</p> 	<p>Compatible with most hard, non-porous surfaces</p>	<ol style="list-style-type: none"> 1. Pre-clean heavily soiled areas. 2. Wipe hard, non-porous environmental surfaces with towelette. 3. Allow surface to remain wet for 1 minute. For certain listed viruses, allow to remain wet for 30 seconds. 4. Allow to air dry or rinse with potable water if necessary. <p><i>Note: All food contact surfaces must be rinsed with potable water. Do not use on glassware, utensils, or dishes</i></p>
<p>Clorox Healthcare® Bleach Germicidal Disinfectant Wipes</p> 	<p>Autoclaves, blood glucose monitors, IV poles, exam tables, bed railings, dialysis machines, respiratory equipment, ultrasound transducers and probes, gurneys, mammography equipment, plastic mattress covers, phlebotomy trays, reception counters, stretchers, wash basins and wheelchairs</p>	<ol style="list-style-type: none"> 1. Put on PPE 2. Remove gross soil if present or if disinfecting for <i>C. difficile</i> spores 3. Wipe the surface until completely wet 4. Wait for the contact time (3 minutes for all pathogens listed on the product label) 5. Discard the wipe
<p>Super Sani-Cloth</p> 	<p>May be used on hard nonporous surfaces of bed railings, blood glucose meters (glucometers), cabinets, carts, chairs, counters, dental unit instruments trays, exam tables, gurneys, isolettes, IV poles, stethoscopes, stretchers, tables, telephones, toilet seats, diagnostic equipment, patient monitoring equipment, patient support and delivery equipment, and many more! Please refer to the product/master label to get the full listing of use sites.</p>	<ol style="list-style-type: none"> 1. Always dispense wipe through lid. When not in use, keep lid closed to prevent moisture loss. 2. If present, remove heavy soil loads prior to disinfecting. 3. Unfold a clean wipe and thoroughly wet surface. 4. Allow treated surface to remain wet for 2 minutes. Let air dry. 5. Do not reuse towelettes. Dispose of used towelette in trash. DO NOT FLUSH IN TOILET.

Chemical	Used On	How
<p>Sani-Cloth AF3</p> 	<p>ALCOHOL FREE. May be used on hard nonporous surfaces of bed railings, blood glucose meters (glucometers), cabinets, carts, chairs, counters, dental unit instruments trays, exam tables, gurneys, isolettes, IV poles, stethoscopes, stretchers, tables, telephones, toilet seats, diagnostic equipment, patient monitoring equipment, patient support and delivery equipment, and many more! Please refer to the product/master label to get the full listing of use sites</p>	<ol style="list-style-type: none"> 1. Always dispense wipe through lid. When not in use, keep lid closed to prevent moisture loss. 2. If present, remove heavy soil loads prior to disinfecting. 3. Unfold a clean wipe and thoroughly wet surface. 4. Allow treated surface to remain wet for 3 minutes. Let air dry. 5. Do not reuse towelettes. Dispose of used towelette in trash. DO NOT FLUSH IN TOILET.
<p>PDI Easy Screen Cleaning Wipes</p> 	<p>Removes dirt and grime from touch screen surfaces.</p>	<ol style="list-style-type: none"> 1. Pull a single wipe from the canister. Close lid to prevent moisture loss 2. Thoroughly wipe screen's surface. 3. Dispose of used wipe in trash. 4. Do not flush in toilet. 5. Consult Equipment Manufacturer's instructions prior to use to ensure product is compatible with device being cleaned.

For any additional questions or concerns regarding approved chemicals for use on COVID-19, please contact Infection Prevention and Control at 419-383-5006