

NURSING SERVICES GUIDELINE HEMODIALYSIS

Guideline: Phoenix Meter Procedure

Policy Number Superseded:

Responsibility: Hemodialysis Registered Nurse (RN)

Purpose of Guideline:

Equipment: Phoenix Meter



Effective Date:
June 2024

Initial Effective Date:
July 2005

Procedure	Point of Emphasis
(1) Daily QC should be performed and documented prior to use.	To ensure accurate results.
(2) Monthly calibration of Phoenix Meter will be performed and documented by Biomed.	To ensure calibration within manufacturer's recommendations.
(3) Use Phoenix Meter only after changing the dialysis bath.	To view conductivity and temperature of machine.
(4) Connect meter to sample port on venous line.	Make sure that there is no air in the syringe.
(5) Aspirate fluid by pulling the syringe plunger to draw fluid through the device. While the solution is flowing, observe the temperature, conductivity and Ph result.	Make sure the solution is flowing when obtaining the temperature result. Idle solution may cause inaccurate results.
(6) Document readings.	Document Phoenix meter readings and machine readings in appropriate space on Hemodialysis Treatment Record Value Ranges: Conductivity: 13.4-14.0 PH: 6.9-7.6 Temperature: 34.5- 38.5 degrees Celsius

Guideline:
Phoenix Meter Procedure

2

	If reading is out of range, an air bubble may be present, Disconnect the meter from the sample port and repeat reading. Document if readings are in range or out of range. If out of range, notify biomed and do not initiate treatment.
--	--

References: Phoenix Machine Operators Manual

Approved by:

*Kurt Kless, MSN, MBA, RN, NE-BC
Chief Nursing Officer*

Initial effective date:

July 2005

Review/Revision Completed by:

*Revised by: Paula Kosmider, BSN, RN
Reviewed by Policy & Standard
Committee
March 2011, January 2017, January
2020, December 2022*

Review/Revision Date:

April 2008
February 2011
March 28, 2011
May 2014
June 2015
May 2018
April 2021
June 2024

Next review date:

June 2027