


Name of Policy:	<u>Continual Quality Assurance for Megavoltage Radiation Treatment Units</u>	
Policy Number:	3364-134-32	
Department:	Radiation Oncology	
Approving Officer:	Chief Executive Officer, UTMC Clinical Service Chief, Radiation Oncology	
Responsible Agent:	Manager, Radiation Oncology	
Scope:	Radiation Oncology	Effective Date: 3/1/2021 Initial Effective Date: 9/1980
<input type="checkbox"/> New policy proposal <input checked="" type="checkbox"/> Minor/technical revision of existing policy <input type="checkbox"/> Major revision of existing policy <input type="checkbox"/> Reaffirmation of existing policy		

(A) Policy Statement

Under the supervision of a Board Certified Medical Physicist, each megavoltage treatment machine shall be checked daily, and monthly for physical and mechanical integrity and performance. For a continuous quality assurance process for each treatment modality records of data acquisition are maintained and trend analysis on quality assurance data is performed and reviewed by a QMP.

(B) Purpose of Policy

To ensure satisfactory performance of megavoltage radiation therapy units.

(C) Procedure

- 1) Under the supervision of a Board Certified Medical Physicist, daily, and monthly QA is performed using generally accepted Quality Assurance procedures, based on the AAPM TG 142 report:
 - a) Daily checks are performed prior to patient's treatment and include the following:
 - i) Dosimetry
 - (1) X-ray output constancy
 - (2) Electron output constancy
 - ii) Imaging
 - (1) Collision interlocks
 - (2) Positioning/Repositioning
 - (3) Imaging/Treatment coordinate coincidence
 - (4) Image quality test
 - iii) Mechanical
 - (1) Laser localization
 - (2) Distance indicator (ODI) @ iso
 - (3) Collimator size indicator
 - (4) MLC - Picket fence test
 - iv) Safety
 - (1) Door interlock beam off
 - (2) Door closing safety
 - (3) Audiovisual monitors
 - (4) Radiation area monitor
 - (5) Beam on indicator
 - v) Any other parameter deemed necessary.
 - b) Monthly checks are performed and include the following:
 - i) Dosimetry
 - (1) X-ray output constancy
 - (2) Electron output constancy
 - (3) Backup monitor chamber constancy

- (4) Typical dose rate output constancy
- (5) Photon beam profile constancy
- (6) Electron beam profile constancy
- (7) Electron beam energy constancy
- (8) Dynamic Wedge -- wedge factor for all energies
- ii) Mechanical
 - (1) Light/radiation field coincidence
 - (2) Light/radiation field coincidence (asymmetric)
 - (3) Distance check device for lasers compared with front pointer
 - (4) Gantry/collimator angle indicators
 - (5) Jaw position indicators (asymmetric)
 - (6) Cross hair centering (walkout)
 - (7) Treatment couch position indicators
 - (8) Localizing lasers
- iii) Safety
 - (1) Laser Guard interlock test
- iv) Respiratory Gating
 - (1) Beam output constancy
 - (2) Phase, amplitude beam control
 - (3) In-room respiratory monitoring system
 - (4) Gating interlock
- v) EDW
 - (1) Wedge factor for all energies (45 or 60 deg)
- vi) MLC
 - (1) Setting vs. Radiation Field for two patterns.
 - (2) Travel speed
 - (3) Leaf positional accuracy at all cardinal gantries
- vii) Imaging
 - (1) Planar MV
 - (a) Imaging and treatment coordinate coincidence (four cardinal angles)
 - (b) Scaling
 - (c) Spatial Resolution
 - (d) Contrast
 - (e) Uniformity
 - (f) Noise
 - (2) Planar kV
 - (a) Imaging and treatment coordinate coincidence (four cardinal angles)
 - (b) Scaling for 8 mm diameter disc
 - (c) Spatial Resolution
 - (d) Contrast
 - (e) Uniformity and Noise
 - (3) Cone-Beam CT (kV and MV)
 - (a) Geometric distortion
 - (b) Spatial Resolution
 - (c) Contrast
 - (d) HU constancy
 - (e) Uniformity and Noise
- viii) Additional in-house tests (not required by TG142)
 - (1) Photon energy constancy
- ix) Any other parameter deemed necessary.

- 2) Document and file the results in common drive.
- 3) Under supervision of board certified medical physicists, annual QA will be performed based on recommendations of TG 142 of the AAPM within 14 months of initial acceptance or the previous annual QA. Reports of annual accelerator QA's are kept in common drive.

<p>Approved by:</p> <p><u>/s/</u> _____ <u>03/09/2021</u> Mersiha Hadziahmetovic MD Clinical Service Chief, Radiation Oncology Date</p> <p><u>/s/</u> _____ <u>03/12/2021</u> Richard P. Swaine Chief Executive Officer -UTMC Date</p> <p><i>Review/Revision Completed By: Michelle Giovanoli</i></p>	<p>Review/Revision Date:</p> <p>10/1981 2/1996 1982 3/1997 1983 4/1999 1984 3/2002 1985 1/2007 1986 4/22/2010 1987 10/1/2013 1988 8/1/2016 1989 10/1/2017 2/1990 2/1/2018 1/1994 2/1/2021 2/1995</p> <hr/> <p>Next Review Date: 2/1/2024</p>
<p>Policies Superseded by This Policy: 38-32</p>	