


<b>Name of Policy:</b> <u>Weight Limits of Imaging Tables</u> <b>Policy Number:</b> 3364-135-001 <b>Department:</b> Radiology <b>Approving Officer:</b> Chief Operating Officer - UTMC <b>Responsible Agent:</b> Chairman & Professor, Radiology <b>Scope:</b> Radiology	  <b>Effective Date:</b> 01/01/2023 Initial Effective Date: 06/01/2017
<input type="checkbox"/> New policy proposal <input type="checkbox"/> Minor/technical revision of existing policy <input type="checkbox"/> Major revision of existing policy <input checked="" type="checkbox"/> Reaffirmation of existing policy	

**(A) Policy Statement**

To ensure no patient is placed on the Imaging Tables that exceeds the manufactures weight limit for safe operation.

**(B) Purpose of Policy**

To comply with the Safe Medical Device Act of 1990 and provide a safe procedural environment for all patients of the Cardiovascular/DP Lab.

**(C) Procedure**

1. The Table weight limits are as follows:
  - a. GE Discovery XR650 elevating table 220kg (485 lbs)
  - b. GE Diffinium 8000 elevating table 200 kg (485 lbs)
  - c. GE RF P180 Footboard and Table with tabletop movement = 585 lbs (266 kg)
    - Table without tabletop movement = 710 lbs (323 kg)
  - d. Philips MD Eleva 185 kg (407 lbs)
  - e. GE Prodigy 350 lbs (dexa)
  - f. Toshiba Aquilion One, 320 300 kg (660 lbs)
  - g. Toshiba Aquilion 64 204 kg (450 lbs) Abdominal exams 227 kg (500) Head exams
  - h. GE Premier 3T – 250 kg (550) lbs
  - i. Siemens Espree 1.5T – 200 kg (440 lbs)
  - j. GE Infinia Hawkeye 182 kg (400 lbs)
  - k. Philips Gemini 64 227 kg (500 lbs)
  
2. Patient weight will be obtained from the patient’s record. If patients’ weight is within 15 kg (30 lbs.) of the table limit:
  - a. Patient will be weighted prior to scheduled procedure for a real time accurate weight.
  - b. If patient found to be under the weight limit of the table, the procedure will continue as scheduled.
  - c. If patient is found to be over the weight limit the procedure physician will be notified and the procedure will be rescheduled or cancelled.

<b>Approved by:</b>	<b>Review/Revision Date:</b>
<u>/s/</u> Ryan Landis Director, Radiology	06/01/2017 6/1/2020 01/01/2023
<u>/s/</u> Haitham Elsamaloty, MD Chairman & Professor, Radiology	<u>01/06/2022</u> Date
<u>/s/</u> Christine Stesney-Ridenour, FACHE Chief Operating Officer – UTMC	<u>01/10/2023</u> Date
<u>/s/</u> Christine Stesney-Ridenour, FACHE Chief Operating Officer – UTMC	<u>01/11/2023</u> Date
	Next Review Date: 01/01/2026