

STERILE PROCESSING PROCEDURE STERILIZATION

Procedure SP3-9: BI test for gas plasma sterilizers (STERRAD)



Policy Number Superseded:

Responsibility: Operations Manager,
Sterile Processing

Effective Date:
December 2025

Initial Effective Date:
1996

Purpose of Procedure: To challenge the efficiency of the sterilizer to attain a maximum level of sterility assurance. To follow the IFU of the STERRAD 100NX.

Procedure:

The sterile processing department will perform a biological test on the STERRAD sterilizers each morning on the standard cycle, and the first cycle of any of the following cycles: express, flex or duo.

- (A) Each box of BI test vials contains 30 vials of test media. One vial is to be used as the control each day for the two sterilizers.
- (B) Each cycle can use the same control from the corresponding lot number.
- (C) Biological testing is to be done for the standard cycle, run after midnight in each sterilizer.
- (D) Biological testing is done with the first cycle of the day for express, flex or duo as needed.
- (E) Remove vials from box. Label 1 vial C-5, 1 vial C-6 on the label for control, with date and cycle number. Press top of vial down, crush and place in incubator for control.
- (F) Prepare test vials labeled S-5, S-6 with date and cycle number and place each vial inside STERRAD peel pouch with an indicator.

Procedure:

BI test for gas plasma sterilizers (STERRAD)

2

- (G) Label F for Flex, E for Express and D for Duo with the cycle number and date.
- (H) Enter all pertinent information in the SPM computer system.
- (I) Place the bio-test pack on the bottom shelf near the back of the sterilizer.
- (J) Run load cycle.
- (K) Upon completion of cycle, remove processed pack. Crush vials and place in the designated incubator.
- (L) If the cycle cancels, the BI test pack should be discarded and a new pack used when starting the next cycle.
- (M) Rewrap and re-sterilize any other items which were processed in the cancelled cycle.
- (N) If a true positive result is observed, repeat the biological test procedure on a second vial that has been processed in the sterilizer. If a positive result is again obtained, have the sterilizer checked by bio-med personnel. Follow the item recall procedure.
- (O) Use one control vial of the same lot number per day.
- (P) Read and record results at 24 hours.
- (Q) Scan the cycle printout and enter biological test information and results in the computerized sterilization management system.

Approved by:
Christopher Lehnert
Operations Supervisor, Sterile Processing

Initial effective date:
1996

Review/Revision Completed by:
Operations Supervisor, Sterile Processing

Review/Revision Date:
1998
1999
2002
2005
2007
2009
2010

Procedure:
BI test for gas plasma sterilizers (STERRAD)

3

January 2014
2017
December 2019
December 2025

Next review date:
December 2028