

Full Tanks

- Full tanks will continue to be delivered to the tank room



DOWL 0140

Partially Full

- According to Joint Commission, as soon as tanks or cylinders leave the tank room they are considered **“Partially Full”**



In Use

- Tanks are considered **“in use”** when on a stretcher, a bed, a wheelchair, or code cart.



Empty



- When a cylinder is “**empty**” (below 500 psi), they must be placed in a carrier.
- Carriers are located on each floor in the room by the elevator.
- 3rd Floor “empty” carrier is located in 3145.

Filled Carrier

- When the carrier is filled with “**Empty**” cylinders, it **MUST** be taken to the tank room on Dowling Hall Dock, Room 0140
- Place “**Empty**” cylinders in appropriate location in tank room



Labeling

- The empty cylinder holder will be labeled “**Empty Cylinders**”
- Partially full cylinder holders will be labeled “**Partially Full Cylinders**”
- **IMPORTANT:** Do not commingle “Empty Cylinders” and “Partially Full Cylinders” on the floors
 - This helps eliminate confusion and delay in care
- Transportation carts should be labeled with a “partially Full” sticker. Do not leave empty cylinders in carts

