

Animal Hazard

Complete the form and submit to **DLAR 5** days prior to any hazardous use.



Protocol Number:		Special Instructions
PI: Agent: Biological Hazard 2 CHEM 2 CHEM 3 Start:		 CHEM 2/BSL2: N-95 Mask* and Safety Goggles are mandatory when handling hazards and/ or soiled equipment outside of the BSC. *N-95 respirator requires fit testing through EHRS CHEM 3: Must be handled in a BSC and housed on a negative flow ventilated rack. (Classification of CHEM 2/3 will return to a CHEM 0 seven days after the last treatment and after cage change. ABSL2 remains the same for the duration of experiment.)
# Of cages:	IBC Protocol Number	Sanitation Guidelines
Lab Contact: DLAR: 419-383x4310	Phone: EHRS: 419-530x3600	 Biohazard Cage Cleaning inside the BSC: Empty hazardous bedding into an autoclavable hazard bag in the BSC. The empty hazardous enclosures is bagged and with Hazardous sticker. Spray outside of bag
Personal Protective Equipment Requirements (working inside the Biosafety Cabinet "BSC") • Hands: Disposable Synthetic Gloves □ If checked 2 nd pair of gloves required when handling the hazard. • Face: Surgical Mask • Body: Disposable Gown (single use) □ If checked, Chemo gowns must be worn during prep/administration of hazard • Head: Disposable Bouffant Cap • Shoes: Disposable Shoe Covers		 with disinfectant before removing from BSC/room and transport to the autoclave and autoclave. <u>Chemical hazard Cage Cleaning inside the BSC</u>: Empty hazardous bedding into an autoclavable hazard bag in the BSC. Seal the bag and place in hazardous bin. The empty hazardous enclosures are sprayed with disinfectant before covering and transporting to cage wash. Equipment, surfaces and supplies: Spray with disinfectant prior to removal from room to go to cage wash. Water Disposal: Pour into bedding to be absorbed. Euthanasia & Carcass Disposal: Euthanized in room, placed in hazard bag, label with hazardous sticker. Seal closed and spray with disinfectant to remove from room. Carcasses are place in hazard box on the left side of the walk in cooler in Rm 91. <u>Room Cleaning</u>: Clean room as outlined in DLAR SOP with appropriate disinfectants.