Name of Policy:	Assisting with sphenoidal electrodes	~	
Policy Number:	3364-138-15	THE UNIVERSITY OF TOLEDO	
Department:	Neurodiagnostic Services		
Approving Officer:	Senior Hospital Administrator		
Responsible Agent:	Director, Pulmonary Services		
Scope:	The University of Toledo Medical Center Neurodiagnostic Services	Effective Date:6/1/2023Initial Effective Date:3/10/2004	
New policy proposal X Minor/technical revision of existing policy Major revision of existing policy Reaffirmation of existing policy			

Policy Statement

Assisting with the placement of sphenoidal electrodes will be performed as outlined in this policy.

Procedure

- 1. The sphenoidal cart should be fully stocked at all times.
- 2. The nursing office should be called to schedule conscious sedation. The conscious sedation nurse should report 30 minutes prior to physician arrival.
- 3. The neurology resident is given the information for ordering medications needed from pharmacy.
 - a. IV of normal saline to be running at 40ml/hour just prior to procedure.
 - b. 2% Lidocaine Jelly to be applied by technologist 20 minutes prior to start of procedure.
 - c. Ativan 1mg to be given 1 hour prior to procedure by staff RN.
 - d. Fentanyl 50mcg IV to be given during procedure.
- 4. The area over the zygoma must be free of hair, shaved or clipped back out of the way.
- 5. The physician will insert electrodes using sterile technique.
- 6. Once inserted and secured, the technologist will apply dressing and attach the connectors.
- 7. The electrodes will be plugged into the headbox.
- 8. A sphenoidal montage will be used for rest of recording.

Approved by:		Review/Revision Date:
		09/13/2004 10/03/2006
<u>/s/</u>	6/22/2023	06/29/2007
Michael Taylor	Date	08/11/2010
Director, Pulmonary Services		08/08/2011 08/01/2014
		06/01/2017
/s/	6/29/2023	05/29/2020 - 06/01/2023
Russell Smith	Date	00/01/2025
Senior Hospital Administrator		
		Next Review Date: 6/1/2026