


Name of Policy: <u>LTME/Video EEG Monitoring Checks</u> Policy Number: 3364-138-20 Department: Neurodiagnostics Approving Officer: Associate VP of Patient Care Services UTMC Responsible Agent: Manager, Neurodiagnostics Scope: The University of Toledo Medical Center Neurodiagnostics	 Effective Date: 6/1/2020 Initial Effective Date: 5/27/2014
<input type="checkbox"/> New policy proposal <input type="checkbox"/> Major revision of existing policy	
<input type="checkbox"/> Minor/technical revision of existing policy <input checked="" type="checkbox"/> Reaffirmation of existing policy	

Purpose of this procedure is to establish guidelines for checking running LTME/Video EEG monitoring patients.

All patients on LTME/Video EEG monitoring will be checked at least twice daily for the following.

1. Impedance
2. EKG patches.
3. Check the connection of imaging friendly electrodes.
4. Check for active connection to the server, if available.
5. Make sure EEG is recording.
6. The camera should be position to appropriately visualize patient.
7. Check for any sedation running.
8. Check that cables and cords are not tangled.

Document on the recording:

1. Impedance check.
2. That cable is correctly connected to the breakout or EMU40 box.
3. Sedatives used and dosage.

Approved by: <u>/s/</u> Cynthia Zapotosky, BSN, RN Manager Neurodiagnostics Date _____ <u>/s/</u> Monecca Smith, MSN, RN Associate Vice President of Patient Care Services Date _____	Review/Revision Date: 5/27/2014 5/29/2020 Next Review Date: 6/1/2023
Policies Superseded by This Policy:	